

PACO42394

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC, IF

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200,
IF:1:50-1:200

Protein Background:

GABA, the major inhibitory neurotransmitter in the vertebrate brain, mediates neuronal inhibition by binding to the GABA/benzodiazepine receptor and opening an integral chloride channel.

Gene ID:

GABRQ

Uniprot

Q9UN88

Synonyms:

Gamma-aminobutyric acid, receptor subunit theta (GABA(A) receptor subunit theta), GABRQ


Immunogen:

Recombinant Human Gamma-aminobutyric acid, receptor subunit theta protein (348-611AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human brain tissue using PACO42394 at dilution of 1:100.



Immunofluorescent analysis of MCF-7 cells using PACO42394 at dilution of 1:100 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L).