

# GABRQ Antibody, Biotin conjugated



PACO42397

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:****Protein Background:**

GABA, the major inhibitory neurotransmitter in the vertebrate brain, mediates neuronal inhibition by binding to the GABA/benzodiazepine receptor and opening an integral chloride channel.

**Gene ID:**

GABRQ

**Uniprot**

Q9UN88

**Synonyms:**

Gamma-aminobutyric acid, receptor subunit theta (GABA(A) receptor subunit theta), GABRQ

**Immunogen:**

Recombinant Human Gamma-aminobutyric acid, receptor subunit theta protein (348-611AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A