

HPCAL4 Antibody, FITC conjugated



PACO42436

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

May be involved in the calcium-dependent regulation of rhodopsin phosphorylation.

Gene ID:

HPCAL4

Uniprot

Q9UM19

Synonyms:

Hippocalcin-like protein 4 (HLP4), HPCAL4

Immunogen:

Recombinant Human Hippocalcin-like protein 4 protein (2-191AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A