

KIAA1257 Antibody



PACO42450

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

KIAA1257

Uniprot

Q9ULG3

Synonyms:

Uncharacterized protein KIAA1257, KIAA1257


Immunogen:

Recombinant Human Uncharacterized protein KIAA1257 protein (1-297AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human breast cancer using PACO42450 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO42450 at dilution of 1:100.