

C19orf53 Antibody, Biotin conjugated



PACO42469

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May have a potential role in hypercalcemia of malignancy.

Gene ID:

C19orf53

Uniprot

Q9UNZ5

Synonyms:

Leydig cell tumor 10 kDa protein homolog, C19orf53

Immunogen:

Recombinant Human Leydig cell tumor 10 kDa protein homolog protein (1-99AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A