

# MCTS1 Antibody, FITC conjugated



PACO42480

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

MCTS1

**Uniprot**

Q9ULC4

**Synonyms:**

Malignant T-cell-amplified sequence 1 (MCT-1) (Multiple copies T-cell malignancies), MCTS1, MCT1

**Immunogen:**

Recombinant Human Malignant T-cell-amplified sequence 1 protein (1-181AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A