

# GDF2 Antibody, Biotin conjugated



PACO42501

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Potent circulating inhibitor of angiogenesis. Could be involved in bone formation. Signals through the type I activin receptor ACVRL1 but not other Alks. Signaling through SMAD1 in endothelial cells requires TGF-beta coreceptor endoglin/ENG.

**Gene ID:**

GDF2

**Uniprot**

Q9UK05

**Synonyms:**

Growth/differentiation factor 2 (GDF-2) (Bone morphogenetic protein 9) (BMP-9), GDF2, BMP9

**Immunogen:**

Recombinant Human Growth/differentiation factor 2 protein (300-429AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A