

# PADI2 Antibody, Biotin conjugated



PACO42605

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

Catalyzes the deimination of arginine residues of proteins.

**Gene ID:**

PADI2

**Uniprot**

Q9Y2J8

**Synonyms:**

Protein-arginine deiminase type-2 (EC 3.5.3.15) (PAD-H19) (Peptidylarginine deiminase II) (Protein-arginine deiminase type II), PADI2, KIAA0994 PAD2 PDI2

**Immunogen:**

Recombinant Human Protein-arginine deiminase type-2 protein (1-665AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A