

MRPS28 Antibody



PACO42606

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

mitochondrial inner membrane, mitochondrial small ribosomal subunit, mitochondrion, poly(A) RNA binding, mitochondrial translation, mitochondrial translational elongation, mitochondrial translational initiation, mitochondrial translational termination, organelle organization.

Gene ID:

MRPS28

Uniprot

Q9Y2Q9

Synonyms:

28S ribosomal protein S28, mitochondrial (MRP-S28) (S28mt) (28S ribosomal protein S35, mitochondrial) (MRP-S35) (S35mt) (Mitochondrial small ribosomal subunit protein bS1m), MRPS28, MRPS35


Immunogen:

Recombinant Human 28S ribosomal protein S28, mitochondrial protein (72-187AA).


Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images



Immunohistochemistry of paraffin-embedded human appendix tissue using PACO42606 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human adrenal gland tissue using PACO42606 at dilution of 1:100.