

SETBP1 Antibody



PACO42754

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

SETBP1

Uniprot

Q9Y6X0

Synonyms:

SET-binding protein (SEB), SETBP1, KIAA0437

Immunogen:

Recombinant Human SET-binding protein (1-242AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A