

# TP53TG5 Antibody



PACO42766

---

## Product Information

**Size:**

50ug

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

May play a significant role in p53/TP53-mediating signaling pathway.

**Gene ID:**

TP53TG5

**Uniprot**

Q9Y2B4

**Synonyms:**

TP53-target gene 5 protein (TP53-inducible gene 5 protein), TP53TG5, C20orf10

**Immunogen:**


Recombinant Human TP53-target gene 5 protein (1-290AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---



Immunohistochemistry of paraffin-embedded human heart tissue using PACO42766 at dilution of 1:100.