

TP53TG5 Antibody, HRP conjugated



PACO42767

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

May play a significant role in p53/TP53-mediating signaling pathway.

Gene ID:

TP53TG5

Uniprot

Q9Y2B4

Synonyms:

TP53-target gene 5 protein (TP53-inducible gene 5 protein), TP53TG5, C20orf10

Immunogen:

Recombinant Human TP53-target gene 5 protein (1-290AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A