

NDUFAF1 Antibody, Biotin conjugated



PACO42785

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:**

Chaperone protein involved in the assembly of the mitochondrial NADH: ubiquinone oxidoreductase complex (complex I).

Gene ID:

NDUFAF1

Uniprot

Q9Y375

Synonyms:

Complex I intermediate-associated protein 30, mitochondrial (NADH dehydrogenase [ubiquinone] 1 alpha subcomplex assembly factor 1), NDUFAF1, CIA30

Immunogen:

Recombinant Human Complex I intermediate-associated protein 30, mitochondrial protein (25-327AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A