

C15orf41 Antibody, FITC conjugated



PACO42820

Product Information

Size:

50ug

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:**Protein Background:****Gene ID:**

C15orf41

Uniprot

Q9Y2V0

Synonyms:

Uncharacterized protein C15orf41 (Protein HH114), C15orf41

Immunogen:

Recombinant Human Uncharacterized protein C15orf41 protein (1-183AA).

Storage:

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

Product Images

N/A

N/A