

# Cd81 Antibody



PACO42890

---

## Product Information

**Size:**

50ug

**Reactivity:**

Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

May play an important role in the regulation of lymphoma cell growth. May be involved in the acrosome reaction.

**Gene ID:**

Cd81

**Uniprot**

P35762

**Synonyms:**

CD81 antigen (26 kDa cell surface protein TAPA-1) (Target of the antiproliferative antibody 1) (CD antigen CD81), Cd81, Tapa1

**Immunogen:**

Recombinant Mouse CD81 antigen protein (116-201AA).

**Storage:**

Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol, 0.01M PBS, PH 7.4

## Product Images

---

N/A

N/A