

PDPK1 Antibody



PACO42903

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:20-1:200

Protein Background:

Gene ID:

PDPK1

Uniprot

O15530

Synonyms:

3-phosphoinositide-dependent protein kinase 1 (hPDK1) (EC 2.7.11.1), PDPK1, PDK1


Immunogen:

Recombinant Human 3-phosphoinositide-dependent protein kinase 1 protein (277-556AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot. All lanes: PDPK1 antibody at 2.5 μ g/ml. Lane 1: HeLa whole cell lysate. Lane 2: A549 whole cell lysate. Secondary: Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 64, 59, 61, 49, 51 kDa. Observed band size: 64 kDa.



Immunohistochemistry of paraffin-embedded human thyroid tissue using PACO42903 at dilution of 1:100.