

ATXN1 Antibody



PACO42955

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Chromatin-binding factor that repress Notch signaling in the absence of Notch intracellular domain by acting as a CBF1 corepressor. Binds to the HEY promoter and might assist, along with NCOR2, RBPJ-mediated repression. Binds RNA in vitro. May be involved in RNA metabolism. The expansion of the polyglutamine tract may alter this function.

Gene ID:

ATXN1

Uniprot

P54253

Synonyms:

Ataxin-1 (Spinocerebellar ataxia type 1 protein), ATXN1, ATX1 SCA1


Immunogen:

Recombinant Human Ataxin-1 protein (586-815AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human gallbladder tissue using PACO42955 at dilution of 1:100.