

XIAP Antibody



PACO42968

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:20-1:200

Protein Background:**Gene ID:**

XIAP

Uniprot

P98170

Synonyms:

E3 ubiquitin-protein ligase XIAP (EC 2.3.2.27) (Baculoviral IAP repeat-containing protein 4) (IAP-like protein) (ILP) (hILP) (Inhibitor of apoptosis protein 3) (IAP-3) (hIAP-3) (hIAP3) (RING-type E3 ubiquitin transferase XIAP) (X-linked inhibitor of apoptosis protein) (X-linked IAP), XIAP, API3 BIRC4 IAP3


Immunogen:

Recombinant Human E3 ubiquitin-protein ligase XIAP protein (1-200AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Western blot. All lanes: XIAP antibody at 3 μ g/ml. Lane 1: HepG2 whole cell lysate. Lane 2: MCF-7 whole cell lysate. Lane 3: PC-3 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 57 kDa. Observed band size: 57 kDa.



Immunohistochemistry of paraffin-embedded human melanoma using PACO42968 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO42968 at dilution of 1:100.