

TP53BP1 Antibody



PACO42978

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Plays a key role in the response to DNA damage. May have a role in checkpoint signaling during mitosis. Enhances TP53-mediated transcriptional activation.

Gene ID:

TP53BP1

Uniprot

Q12888

Synonyms:

TP53-binding protein 1 (53BP1) (p53-binding protein 1) (p53BP1), TP53BP1


Immunogen:

Recombinant Human Tumor suppressor p53-binding protein 1 protein (1758-1972AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human colon tissue using PACO42978 at dilution of 1:100.