

CD48 Antibody



PACO43072

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Ligand for CD2. Might facilitate interaction between activated lymphocytes. Probably involved in regulating T-cell activation.

Gene ID:

CD48

Uniprot

P09326

Synonyms:

CD48 antigen (B-lymphocyte activation marker BLAST-1) (BCM1 surface antigen) (Leukocyte antigen MEM-102) (SLAM family member 2) (SLAMF2) (Signaling lymphocytic activation molecule 2) (TCT.1) (CD antigen CD48), CD48, BCM1 BLAST1


Immunogen:

Recombinant Human CD48 antigen protein (27-125AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO43072 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human epityphlon tissue using PACO43072 at dilution of 1:100.