

GOLM1 Antibody



PACO43093

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Unknown. Cellular response protein to viral infection.

Gene ID:

GOLM1

Uniprot

Q8NBJ4

Synonyms:

Golgi membrane protein 1 (Golgi membrane protein GP73) (Golgi phosphoprotein 2), GOLM1, C9orf155 GOLPH2


Immunogen:

Recombinant Human Golgi membrane protein 1 protein (112-401AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human colon tissue using PACO43093 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO43093 at dilution of 1:100.