

AMY1A Antibody



PACO43125

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:**Gene ID:**

AMY1A

Uniprot

P04745

Synonyms:

Alpha-amylase 1 (EC 3.2.1.1) (1,4-alpha-D-glucan glucohydrolase 1) (Salivary alpha-amylase), AMY1A; AMY1B; AMY1C, AMY1; AMY1; AMY1


Immunogen:

Recombinant Human Alpha-amylase 1 protein (16-200AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human salivary gland tissue using PACO43125 at dilution of 1:100.