

CD244 Antibody



PACO43130

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Modulate other receptor-ligand interactions to enhance leukocyte activation. CD244/2B4 is the only heterophilic receptor of SLAM family.

Gene ID:

CD244

Uniprot

Q9BZW8

Synonyms:

Natural killer cell receptor 2B4 (NK cell activation-inducing ligand) (NAIL) (NK cell type I receptor protein 2B4) (NKR2B4) (h2B4) (SLAM family member 4) (SLAMF4) (Signaling lymphocytic activation molecule 4) (CD antigen CD244), CD244, 2B4


Immunogen:

Recombinant Human Natural killer cell receptor 2B4 protein (22-224AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human small intestine tissue using PACO43130 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO43130 at dilution of 1:100.