

# DYRK2 Antibody



PACO43181

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

**Gene ID:**

DYRK2

**Uniprot**

Q92630

**Synonyms:**

Dual specificity tyrosine-phosphorylation-regulated kinase 2 (EC 2.7.12.1), DYRK2

**Immunogen:**


Recombinant Human Dual specificity tyrosine-phosphorylation-regulated kinase 2 protein (502-601AA).

**Storage:**


PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Immunohistochemistry of paraffin-embedded human breast cancer using PACO43181 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO43181 at dilution of 1:100.