

PACO43184

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Catalyzes a salvage reaction resulting in the formation of AMP, that is energetically less costly than de novo synthesis.

**Gene ID:**

APRT

**Uniprot**

P07741

**Synonyms:**

Adenine phosphoribosyltransferase (APRT) (EC 2.4.2.7), APRT

**Immunogen:**


Recombinant Human Adenine phosphoribosyltransferase protein (1-180AA).

**Storage:**

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Immunohistochemistry of paraffin-embedded human liver cancer using PACO43184 at dilution of 1:100.