

DPP8 Antibody



PACO43192

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Dipeptidyl peptidase that cleaves off N-terminal dipeptides from proteins having a Pro or Ala residue at position 2. May play a role in T-cell activation and immune function.

Gene ID:

DPP8

Uniprot

Q6V1X1

Synonyms:

Dipeptidyl peptidase 8 (DP8) (EC 3.4.14.5) (Dipeptidyl peptidase IV-related protein 1) (DPRP-1) (Dipeptidyl peptidase VIII) (DPP VIII) (Prolyl dipeptidase DPP8), DPP8, DPRP1


Immunogen:

Recombinant Human Dipeptidyl peptidase 8 protein (10-230AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human endometrial cancer using PACO43192 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human testis tissue using PACO43192 at dilution of 1:100.