

IGFBP4 Antibody



PACO43201

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

IGF-binding proteins prolong the half-life of the IGFs and have been shown to either inhibit or stimulate the growth promoting effects of the IGFs on cell culture. They alter the interaction of IGFs with their cell surface receptors.

Gene ID:

IGFBP4

Uniprot

P22692

Synonyms:

Insulin-like growth factor-binding protein 4 (IBP-4) (IGF-binding protein 4) (IGFBP-4), IGFBP4, IBP4


Immunogen:

Recombinant Human Insulin-like growth factor-binding protein 4 protein (22-258AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human colon cancer using PACO43201 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO43201 at dilution of 1:100.