

PACO43210

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000

Protein Background:

Gene ID:

PIAS2

Uniprot

O75928

Synonyms:

E3 SUMO-protein ligase PIAS2 (EC 2.3.2. -) (Androgen receptor-interacting protein 3) (ARIP3) (DAB2-interacting protein) (DIP) (E3 SUMO-protein transferase PIAS2) (Msx-interacting zinc finger protein) (Miz1) (PIAS-NY protein) (Protein inhibitor of activated STAT x) (Protein inhibitor of activated STAT2), PIAS2, PIASX


Immunogen:

Recombinant Human E3 SUMO-protein ligase PIAS2 protein (532-621AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot. All lanes: PIAS2 antibody at 2 μ g/ml. Lane 1: A431 whole cell lysate. Lane 2: HepG2 whole cell lysate. Lane 3: Mouse brain tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 69, 64, 46 kDa. Observed band size: 69 kDa.