

ADAMTSL1 Antibody



PACO43269

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:**Gene ID:**

ADAMTSL1

Uniprot

Q8N6G6

Synonyms:

ADAMTS-like protein 1 (ADAMTSL-1) (Punctin-1), ADAMTSL1, ADAMTSR1 C9orf94


Immunogen:

Recombinant Human ADAMTS-like protein 1 protein (1483-1762AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human testis tissue using PACO43269 at dilution of 1:100.