

# IRAK4 Antibody



PACO43280

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

IRAK4

**Uniprot**

Q9NWZ3

**Synonyms:**

Interleukin-1 receptor-associated kinase 4 (IRAK-4) (EC 2.7.11.1) (Renal carcinoma antigen NY-REN-64), IRAK4

**Immunogen:**

Recombinant Human Interleukin-1 receptor-associated kinase 4 protein (1-220AA).

**Storage:**

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---

N/A

N/A