

HMGN2 Antibody



PACO43288

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Binds to the inner side of the nucleosomal DNA thus altering the interaction between the DNA and the histone octamer. May be involved in the process which maintains transcribable genes in a unique chromatin conformation.

Gene ID:

HMGN2

Uniprot

P05204

Synonyms:

Non-histone chromosomal protein HMG-17 (High mobility group nucleosome-binding domain-containing protein 2), HMGN2, HMG17


Immunogen:

Recombinant Human Non-histone chromosomal protein HMG-17 protein (1-90AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human thyroid tissue using PACO43288 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human colon cancer using PACO43288 at dilution of 1:100.