

MAPKAP1 Antibody



PACO43291

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000

Protein Background:**Gene ID:**

MAPKAP1

Uniprot

Q9BPZ7

Synonyms:

Target of rapamycin complex 2 subunit MAPKAP1 (TORC2 subunit MAPKAP1) (Mitogen-activated protein kinase 2-associated protein 1) (Stress-activated map kinase-interacting protein 1) (SAPK-interacting protein 1) (mSIN1), MAPKAP1, MIP1 SIN1


Immunogen:

Recombinant Human Target of rapamycin complex 2 subunit MAPKAP1 protein (1-290AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot. All lanes: Target of rapamycin complex 2 subunit
MAPKAP1 antibody at 2 μ g/ml. Lane 1: Hela whole cell lysate. Lane 2:
MCF-7 whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at
1/10000 dilution. Predicted band size: 60, 55, 54, 38, 37, 42 kDa.
Observed band size: 60 kDa.