

PAX7 Antibody



PACO43298

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:200-1:1000,
IHC:1:20-1:200

Protein Background:

Transcription factor playing a role in myogenesis through regulation of muscle precursor cells proliferation.

Gene ID:

PAX7

Uniprot

P23759

Synonyms:

Paired box protein Pax-7 (HuP1), PAX7, HUP1


Immunogen:

Recombinant Human Paired box protein Pax-7 protein (276-505AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot. All lanes: PAX7 antibody at 4 μ g/ml + 293T whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 56, 57 kDa. Observed band size: 56 kDa.



Immunohistochemistry of paraffin-embedded human lung cancer using PACO43298 at dilution of 1:100.