

SOX14 Antibody



PACO43321

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Acts as a negative regulator of transcription.

Gene ID:

SOX14

Uniprot

O95416

Synonyms:

Transcription factor SOX-14 (Protein SOX-28), SOX14, SOX28


Immunogen:

Recombinant Human Transcription factor SOX-14 protein (91-240AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human glioma using PACO43321 at dilution of 1:100.