
Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Mediates in vitro the transfer of all common phospholipids, cholesterol and gangliosides between membranes. May play a role in regulating steroidogenesis.

Gene ID:

SCP2

Uniprot

P22307

Synonyms:

Non-specific lipid-transfer protein (NSL-TP) (EC 2.3.1.176) (Propanoyl-CoA C-acyltransferase) (SCP-chi) (SCPX) (Sterol carrier protein 2) (SCP-2) (Sterol carrier protein X) (SCP-X), SCP2


Immunogen:

Recombinant Human Non-specific lipid-transfer protein (1-322AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO43331 at dilution of 1:100.