

HSPB8 Antibody



PACO43334

Product Information

Size:

50ul

Reactivity:

Human, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000

Protein Background:

Displays temperature-dependent chaperone activity.

Gene ID:

HSPB8

Uniprot

Q9UJY1

Synonyms:

Heat shock protein beta-8 (HspB8) (Alpha-crystallin C chain) (E2-induced gene 1 protein) (Protein kinase H11) (Small stress protein-like protein HSP22), HSPB8, CRYAC E2IG1 HSP22


Immunogen:

Recombinant Human Heat shock protein beta-8 protein (1-196AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot. All lanes: Heat shock protein beta-8 antibody at 4 μ g/ml + Rat heart tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 22 kDa. Observed band size: 22 kDa.