

PACO43342

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

May play a role in the regulation of plasmin-mediated matrix remodeling. Inhibits trypsin, plasmin, factor VIIa/tissue factor and weakly factor Xa. Has no effect on thrombin.

Gene ID:

TFPI2

Uniprot

P48307

Synonyms:

Tissue factor pathway inhibitor 2 (TFPI-2) (Placental protein 5) (PP5), TFPI2


Immunogen:

Recombinant Human Tissue factor pathway inhibitor 2 protein (23-235AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Immunohistochemistry of paraffin-embedded human breast cancer using PACO43342 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO43342 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human tonsil tissue using PACO43342 at dilution of 1:100.