

CTNNBIP1 Antibody



PACO43362

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA

Recommended dilutions:

Protein Background:

Prevents the interaction between CTNNB1 and TCF family members, and acts as negative regulator of the Wnt signaling pathway.

Gene ID:

CTNNBIP1

Uniprot

Q9NSA3

Synonyms:

Beta-catenin-interacting protein 1 (Inhibitor of beta-catenin and Tcf-4), CTNNBIP1, ICAT

Immunogen:

Recombinant Human Beta-catenin-interacting protein 1 protein (1-81AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images

N/A

N/A