

# AMPD3 Antibody



PACO43427

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

AMP deaminase plays a critical role in energy metabolism.

**Gene ID:**

AMPD3

**Uniprot**

Q01432

**Synonyms:**

AMP deaminase 3 (EC 3.5.4.6) (AMP deaminase isoform E) (Erythrocyte AMP deaminase), AMPD3

**Immunogen:**


Recombinant Human AMP deaminase 3 protein (1-260AA).

**Storage:**


PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO43427 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human spleen tissue using PACO43427 at dilution of 1:100.