
Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Converts lysophosphatidic acid, (LPA) into phosphatidic acid, by incorporating an acyl moiety at the sn-2 position of the glycerol backbone.

Gene ID:

AGPAT1

Uniprot

Q99943

Synonyms:

1-acyl-sn-glycerol-3-phosphate acyltransferase alpha (EC 2.3.1.51) (1-acylglycerol-3-phosphate O-acyltransferase 1) (1-AGP acyltransferase 1) (1-AGPAT 1) (Lysophosphatidic acid, acyltransferase alpha) (LPAAT-alpha) (Protein G15), AGPAT1, G15


Immunogen:

Recombinant Human 1-acyl-sn-glycerol-3-phosphate acyltransferase α protein (204-283AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO43438 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO43438 at dilution of 1:100.