

BAP1 Antibody



PACO43441

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:200-1:1000

Protein Background:**Gene ID:**

BAP1

Uniprot

Q92560

Synonyms:

Ubiquitin carboxyl-terminal hydrolase BAP1 (EC 3.4.19.12) (BRCA1-associated protein 1) (Cerebral protein 6), BAP1, KIAA0272


Immunogen:

Recombinant Human Ubiquitin carboxyl-terminal hydrolase BAP1 protein (460-729AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot. All lanes: BAP1 antibody at 8 μ g/ml + Hela whole cell lysate. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 80 kDa. Observed band size: 80 kDa.