

RNF125 Antibody



PACO43476

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

E3 ubiquitin-protein ligase that acts as a positive regulator of T-cell activation. E3 ligase proteins mediate ubiquitination and subsequent proteasomal degradation of target proteins.

Gene ID:

RNF125

Uniprot

Q96EQ8

Synonyms:

E3 ubiquitin-protein ligase RNF125 (EC 2.3.2.27) (RING finger protein 125) (T-cell RING activation protein 1) (TRAC-1), RNF125


Immunogen:

Recombinant Human E3 ubiquitin-protein ligase RNF125 protein (1-232AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human skin tissue using PACO43476 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human spleen tissue using PACO43476 at dilution of 1:100.