

GSTA2 Antibody



PACO43521

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000,
IHC:1:20-1:200

Protein Background:

Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles.

Gene ID:

GSTA2

Uniprot

P09210

Synonyms:

Glutathione S-transferase A2 (EC 2.5.1.18) (GST HA subunit 2) (GST class-alpha member 2) (GST-gamma) (GSTA2-2) (GTH2), GSTA2, GST2


Immunogen:

Recombinant Human Glutathione S-transferase A2 protein (1-222AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.


Product Images



Western blot. All lanes: Glutathione S-transferase A2 antibody at 4 μ g/ml + Mouse liver tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 26 kDa. Observed band size: 26 kDa.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO43521 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human liver cancer using PACO43521 at dilution of 1:100.