

MAPK13 Antibody



PACO43524

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Gene ID:

MAPK13

Uniprot

O15264

Synonyms:

Mitogen-activated protein kinase 13 (MAP kinase 13) (MAPK 13) (EC 2.7.11.24) (Mitogen-activated protein kinase p38 delta) (MAP kinase p38 delta) (Stress-activated protein kinase 4), MAPK13, PRKM13 SAPK4


Immunogen:

Recombinant Human Mitogen-activated protein kinase 13 protein (1-365AA).


Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human prostate tissue using PACO43524 at dilution of 1:100.



Immunohistochemistry of paraffin-embedded human kidney tissue using PACO43524 at dilution of 1:100.