

# HOXB1 Antibody



PACO43543

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

ELISA:1:2000-1:10000, WB:1:200-1:1000

**Protein Background:**

Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis. Acts on the anterior body structures.

**Gene ID:**

HOXB1

**Uniprot**

P14653

**Synonyms:**

Homeobox protein Hox-B1 (Homeobox protein Hox-2I), HOXB1, HOX2I

**Immunogen:**


Recombinant Human Homeobox protein Hox-B1 protein (35-180AA).

**Storage:**

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Western blot. All lanes: Homeobox protein Hox-B1 antibody at 9 $\mu$ g/ml + Mouse liver tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 33, 25 kDa. Observed band size: 33 kDa.