

# PRND Antibody



PACO43550

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

**Gene ID:**

PRND

**Uniprot**

Q9UKY0

**Synonyms:**

Prion-like protein doppel (PrPLP) (Prion protein 2), PRND, DPL

**Immunogen:**


Recombinant Human Prion-like protein doppel protein (20-160AA).

**Storage:**

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Immunohistochemistry of paraffin-embedded human testis tissue using PACO43550 at dilution of 1:100.