

PACO43569

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Inhibits the cytotoxic effects of TNF-alpha and overexpressed TNF receptor adapters TRADD, FADD, and RIPK1. Involved in TGF-beta1 inhibition of I kappa B-alpha expression and suppression of TNF-mediated I kappa B-alpha degradation.

**Gene ID:**

TIAF1

**Uniprot**

O95411

**Synonyms:**

TGFB1-induced anti-apoptotic factor 1 (12 kDa TGF-beta-1-induced antiapoptotic factor), TIAF1

**Immunogen:**


Recombinant Human TGFB1-induced anti-apoptotic factor 1 protein (1-115AA).

**Storage:**

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Immunohistochemistry of paraffin-embedded human skin tissue using PACO43569 at dilution of 1:100.