

# SOAT1 Antibody



PACO43586

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, IHC

**Recommended dilutions:**

ELISA:1:2000-1:10000, IHC:1:20-1:200

**Protein Background:**

Catalyzes the formation of fatty acid, cholesterol esters, which are less soluble in membranes than cholesterol. Plays a role in lipoprotein assembly and dietary cholesterol absorption. In addition to its acyltransferase activity, it may act as a ligase.

**Gene ID:**

SOAT1

**Uniprot**

P35610

**Synonyms:**

Sterol O-acyltransferase 1 (EC 2.3.1.26) (Acyl-coenzyme A: cholesterol acyltransferase 1) (ACAT-1) (Cholesterol acyltransferase 1), SOAT1, AACT ACAT1 ACAT ACAT1 SOAT STAT

**Immunogen:**


Recombinant Human Sterol O-acyltransferase 1 protein (1-130AA).

**Storage:**

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Product Images

---



Immunohistochemistry of paraffin-embedded human placenta tissue using PACO43586 at dilution of 1:100.