

CXXC1 Antibody



PACO43587

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

ELISA:1:2000-1:10000, IHC:1:20-1:200

Protein Background:

Transcriptional activator that exhibits a unique DNA binding specificity for CpG unmethylated motifs with a preference for CpGG.

Gene ID:

CXXC1

Uniprot

Q9P0U4

Synonyms:

CXXC-type zinc finger protein 1 (CpG-binding protein) (PHD finger and CXXC domain-containing protein 1), CXXC1, CFP1, CGBP, PCCX1, PHF18


Immunogen:

Recombinant Human CXXC-type zinc finger protein 1 protein (1-210AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Immunohistochemistry of paraffin-embedded human colon tissue using PACO43587 at dilution of 1:100.