

ERGIC2 Antibody



PACO43607

Product Information

Size:

50ul

Reactivity:

Human, Mouse

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

ELISA:1:2000-1:10000, WB:1:500-1:2000

Protein Background:

Possible role in transport between endoplasmic reticulum and Golgi.

Gene ID:

ERGIC2

Uniprot

Q96RQ1

Synonyms:

Endoplasmic reticulum-Golgi intermediate compartment protein 2, ERGIC2, ERV41
PTX1


Immunogen:

Recombinant Human Endoplasmic reticulum-Golgi intermediate compartment protein 2 protein (55-310AA).

Storage:

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Product Images



Western blot. All lanes: Endoplasmic reticulum-Golgi intermediate compartment protein 2 antibody at 4 μ g/ml. Lane 1: HepG2 whole cell lysate. Lane 2: Jurkat whole cell lysate. Lane 3: Mouse lung tissue. Secondary. Goat polyclonal to rabbit IgG at 1/10000 dilution. Predicted band size: 43 kDa. Observed band size: 43 kDa.